



P:\SW-Tek\Timo en Cleo\Timo 6034\

Signed	KA	Date	02-03-2016
Revised	01 /KA	Date	29-10-2018
Revised	02 /KA	Date	31-03-2020

# Overview parts Timo

Copyright reserved by law      Unit of measurement: mm

Lopital Nederland bv

Type nr: 54006034

## Timo 54006034

Part nr.	QTY.	Art.nr.	Part
1	1	70203501	Frame
2	1	70203509	Backrest
3	1	69608024	Footrest
4	4	70519928	Castor Ø 125 mm braked
5	1	70602007	Seat
6	1	70565001	Armrest
7	4		Snaploc
8	2	75005640	Cylinder head screw M8x16 Super Duplex 4529
9	4		Cap castor
10	4		Insert Ø25
11	4		Countersunk screw M6x10
12	4		Capring M4
13	4		Chipboard screw 4,2x30
14	4		Screw cap M4



P:\SW-Tek\Timo en Cleo\Timo 6034\

Signed	KA	Date	02-03-2016
Revised	01 /KA	Date	29-10-2018
Revised	02 /KA	Date	31-03-2020

Copyright reserved  
by law

Unit of measure-  
ment: mm

## Overview parts Timo

Lopital Nederland bv

Type nr: 54006034